



## 13012 Icarus

---

Multi level stretcher with steel tank



---

### FEATURES

---

Multi level stretcher Icarus is made with stainless steel, has 4 wheels 2 of which with brakes, side rails and foot rest. Four gas pistons allow a perfect adjustment to all necessary heights for a practical and professional work.

Available in following versions with:

- ABS lying part for a perfect cleaning and disinfection (art. 13010)
- stainless steel tank (art. 13012)

---

## MAIN FEATURES

---

|                       |   |
|-----------------------|---|
| Multilevel            |  |
| Lying surface         | Without Trendelenburg   |
| Wheel type            | 2 Rotating Wheels + 2 Small Fixed Wheels  |
| Min height (cm)       | 51  |
| Max height (cm)       | 96  |
| Width (cm)            | 62  |
| Length (cm)           | 198   |
| Weight (kg)           | 39  |
| Loading capacity (kg) | 170   |
| Size                  | Adults  |
| Usage                 | Transport   |
| Material              | Stainless steel   |

---

## SEE ALSO

---



**Icarus T3010**

Multi level stretcher with side rails

---

---

## ACCESSORIES

---



**607/B**

Blue belt in 2 pieces with plastic buckle

---