




## 9050 Nilix

Yellow basic spinal board with pins



### FEATURES

Basic spinal board made of high density plastic material, equipped with 8 PINS. It's supplied with 3 belts.

### MAIN FEATURES

Length (cm)	183
Width (cm)	40
Thickness (cm)	5
Loading capacity (kg)	150
Weight (kg)	6,5
Trauma	
Usage	Immobilization
Material	HDPE (High-density polyethylene)

---

## ACCESSORIES

---



**Jelly 629**

Orange standard head immobilizer  
for spinal board

---