



## 11128/FA

---

On board medical suction unit with autoclavable canister



---

### FEATURES

- Rescue vehicles On board medical suction unit
- Maximum aspiration flow 30lt/min.
- Maximum aspiration Power 800 mbar.
- 1000 ml Autoclavabile canister with disposable (model LN)
- 1000 ml autoclavable canister with safety valve and protection filter directly fixed onto the cap (model FA)
- Dimension cm. 9x8,5x6 depth.
- In compliance with all existing norms.

---

### SEE ALSO

---



**11128/LN**

On board medical suction unit with  
disposable liner

---