



7080/4RG Proof Mercury Lite

Ambulance cot in aluminium with variable height and four swivel castors



FEATURES

Our "Mercury Lite" stretcher is available in the 4 swivel castors version, ensuring excellent manoeuvrability. This device is ideal for narrow spaces. The two front castors can be unblocked by means of an intuitive mechanism. They can revert to their original non-rotating position. The fastening of the stretcher inside of an emergency vehicle is ensured by our certified product 955 Proof (see "Fixing Systems"). The stretcher is supplied with mattress, straps and object holder.

Device certified according to the European harmonized safety standards UNI EN 1865 – UNI EN 1789.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	<u>✓</u>
Material	Painted aluminium HDPE (High-density polyethylene)
Multilevel	₩
Lying surface	With Trendelenburg
Intermediate height (cm)	72
Max height (cm)	95
Max loading height (cm)	66
Min loading height (cm)	64
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40
Loading capacity (kg)	250
Type	Barella autocaricante
Wheel type	4 Big Rotating Wheels
Usage	Transport
Size	Adults

IMAGES









SEE ALSO



Mercury 7070/4RG PROOF Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury 7072/4RG PROOF Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version

ACCESSORIES



12080

"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001

Monting kit for bariatric platform on self loading stretcher



12128

Adjustable system for headrest of Meber-X and Mercury stretchers



12154

"BEAR PODS" Bariatric containment system for selfloading stretchers



12018N

Telescopic foldable phleboclysis stand for self loading stretchers



955 PROOF

EN 1789 fixing system for Proof roll-in stretchers



12070/BG

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



7013/A

3 pcs. Skay fire-proof electrowelded mattress



7014/A

4 pcs. Skay fire-proof electrowelded mattress



12100/N

Black electrowelded shaped cushion with MeBer logo for electowelded mattresses



12106Mattress for AIR-COM stretchers -Anti-bedsore mattress with alternate segments



939 Instrument support for stretchers



938

Oximeb oxygen cylinder holder



12054

Oximeb Maxi - Oxygen cylinder holder



12084-001

Cover for bariatric wings set



12170

Black mattress for right bariatric wing for self-loading stretchers



12172

Black mattress for left bariatric wing for self-loading stretchers



707Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



hook and MeBer band (L. 155 cm)







12025

Black multi usage object holder for automatic loading stretcher

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12072/BG "COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



12072/BR "COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



12200
Bariatric wings set for Mercury selfloading stretchers complete with mattresses, belts with reel and assembly kit

STANDARD ACCESSORIES



7013/N 3 pcs. Skay fire-proof electrowelded mattress



2 pieces black belt with metalic hook and MeBer band (L. 155 cm)



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12025 Black multi usage object holder for automatic loading stretcher