



# 7082/4RG Proof Mercury Lite

Certified aluminium variable height stretcher high version with 4 rotating wheels



## FEATURES



4 rotating wheels version is available for "Mercury Lite" which grants extremely maneuverability. This product is ideal for narrow spaces. Operating an easy command the two frontal wheels can be unlocked or reset to their original configuration. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "[Fixing Systems](#)"). Supplied complete with mattress and belts. This version has the same features of model 7080/4RG Proof, but it allows to load on ambulances with higher loading height.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

\*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	
Material	Painted aluminium HDPE (High-density polyethylene)
Multilevel	
Lying surface	With Trendelenburg
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40,5
Loading capacity (kg)	250

IMAGES



SEE ALSO



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with  
variable height and four swivel  
castors



**Mercury 7070/4RG PROOF**  
Stainless steel variable height  
stretcher certified with 4 rotating  
wheels



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version

## ACCESSORIES



**12072/BG**  
"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12106**  
Mattress for AIR-COM stretchers -  
Anti-bedsore mattress with  
alternate segments



**707**  
Red 6-handle fireproof cover with  
Me.Ber tape for self-loading  
stretcher



**938**  
Oximeb oxygen cylinder holder



**12054**  
Oximeb Maxi - Oxygen cylinder  
holder



**12084-001**  
Cover for bariatric wings set



**12170**  
Black mattress for right bariatric  
wing for self-loading stretchers



**12172**  
Black mattress for left bariatric  
wing for self-loading stretchers



**7013/A**  
3 pcs. Skay fire-proof electrowelded  
mattress



**7014/A**  
4 pcs. Skay fire-proof electrowelded  
mattress



**12100/N**  
Black electrowelded shaped  
cushion with MeBer logo for  
electrowelded mattresses



**955 PROOF**  
EN 1789 fixing system for Proof roll-  
in stretchers



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs



**12025**  
Black multi usage object holder for  
automatic loading stretcher



**12080**  
"MPB" wide bariatric platform for  
self-loading stretchers Frog  
Mercury and Gecko complete of  
belts with reel and fixing for  
stretcher



**12080-001**  
Monting kit for bariatric platform  
on self loading stretcher



**12154**  
"BEAR PODS" Bariatric  
containment system for self-  
loading stretchers



**12018N**

Telescopic foldable phlebotomy stand for self-loading stretchers



**939**

Instrument support for stretchers



**12128**

Adjustable system for headrest of Meber-X and Mercury stretchers



**12070/BG**

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles




**12070/BR**

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BR**

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**12200**

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

## STANDARD ACCESSORIES



**12025**

Black multi usage object holder for automatic loading stretcher



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self-loading stretchers and sedan chairs



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)