



7082/4RG Proof Mercury Lite

Certified aluminium variable height stretcher high version with 4 rotating wheels



FEATURES

4 rotating wheels version is available for "Mercury Lite" which grants extremely maneuverability. This product is ideal for narrow spaces. Operating an easy command the two frontal wheels can be unlocked or reset to thier original configuration. Fixing of stretcher inside the rescue vehicle is granted by certified system Art. 955 Proof (see "Fixing Systems"). Supplied complete with mattress and belts. This version has the same features of model 7080/4RG Proof, but it allows to load on ambulances with higher loading height.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

10G Certification	₩
Material	Painted aluminium HDPE (High-density polyethylene)
Multilevel	₩
Lying surface	With Trendelenburg
Intermediate height (cm)	75
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	40,5
Loading capacity (kg)	250

IMAGES







SEE ALSO



Mercury Lite 7080/4RG Proof Ambulance cot in aluminium with variable height and four swivel castors



Mercury 7070/4RG PROOF Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version

ACCESSORIES



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12106 Mattress for AIR-COM stretchers -Anti-bedsore mattress with alternate segments



707Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



938 Oximeb oxygen cylinder holder



12054 Oximeb Maxi - Oxygen cylinder holder



12084-001 Cover for bariatric wings set



12170Black mattress for right bariatric wing for self-loading stretchers



12172Black mattress for left bariatric wing for self-loading stretchers



7013/A 3 pcs. Skay fire-proof electrowelded mattress



7014/A 4 pcs. Skay fire-proof electrowelded mattress



12100/N Black electrowelded shaped cushion with MeBer logo for electowelded mattresses



955 PROOF EN 1789 fixing system for Proof rollin stretchers



608/A-MEB2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



2 pieces black belt with metalic hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



12025 Black multi usage object holder for automatic loading stretcher



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001Monting kit for bariatric platform on self loading stretcher



12154
"BEAR PODS" Bariatric
containment system for selfloading stretchers



12018NTelescopic foldable phleboclysis stand for self loading stretchers



939
Instrument support for stretchers



12128Adjustable system for headrest of Meber-X and Mercury stretchers



12070/BG "COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BR "COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



12200
Bariatric wings set for Mercury selfloading stretchers complete with mattresses, belts with reel and assembly kit

STANDARD ACCESSORIES



12025Black multi usage object holder for automatic loading stretcher



7013/N 3 pcs. Skay fire-proof electrowelded mattress



612/N-2/MEBAdjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/N-MEB 2 pieces black belt with metalic hook and MeBer band (L. 155 cm)