



## 9214 Gemini

Cervical collar set 4 pieces with carrying case

Art. 9204 - 9206 - 9208 - 9210



### FEATURES

GEMINI cervical collar is a highly professional instrument indispensable for the emergency operator. Composed by 2 parts, it is positioned in the front and rear part of the neck allowing thus a perfect immobilization of the rachis. Central hole for tracheotomy, easy to be used and adjust thanks to the rip-off fastening. Made with stamped polyethylene coloured to identify various size and E.V.A. washable inner lining fixed by plastic rivets. Available in 6 different size.

Art. 9204 - 9206 - 9208 - 9210

### MAIN FEATURES

Trauma



### SEE ALSO



**Gemini 9204**

Rigid cervical collar in 2 pieces size  
NECKLESS



**Gemini 9206**

Rigid cervical collar in 2 pieces size  
SHORT



**Gemini 9208**

Rigid cervical collar in 2 pieces size  
REGULAR



**Gemini 9210**

Rigid cervical collar in 2 pieces size  
TALL