



## 630/N Maxima

Black painted anti trauma scoop stretcher



### FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

- Effective hooking and unhooking device to open/close the stretcher.
- Adjustable length with automatic locking device.
- Bag with n. 3 quick unhooking safety belts to immobilize the patient (603/MEB).
- Easy arrangement inside the ambulance when folded.
- Completely built with light-weight and sturdy aluminium.

### MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Épaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	
Material	Painted aluminium

---

## SEE ALSO

---



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher



**Maxima 630/B**  
Blue painted anti trauma scoop  
stretcher



**Maxima 630/G**  
Yellow painted anti trauma scoop  
stretcher



**Maxima 630/N**  
Black painted anti trauma scoop  
stretcher



**Maxima 630/N**  
Black painted anti trauma scoop  
stretcher



**Maxima 630/V**  
Green painted anti trauma scoop  
stretcher



**Maxima 631/G**  
Yellow painted anti trauma scoop  
stretcher with head/trunk in one  
piece



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc

---

## ACCESSORIES

---



**616**  
Certified vertical fastening system  
for scoop stretchers



**616/A**  
Certified vertical fixing system for  
scoop stretchers



**614**  
Certified horizontal fastening  
system for scoop stretchers



**630/CS**  
Carrying case for scoop stretchers



**Double 624**  
Head immobilizer for scoop  
stretcher



**603/MEB**  
Nylon bag containing 3 belts Art.  
607/Ameb