



## 9128 - Me.Ber. Vacuum Baby Fix

---

Vacuum mattress for fixing system with belts for transport of children inside the ambulance





---

### FEATURES

Me.Ber. Vacuum Baby Fix is a vacuum mattress to be used with Art. 9130. Extreme lightness and sturdiness, security belts and perfect encircling of patient make Art. 9128 unique. **Certified UNI EN 1789**, X-ray transparent, nonmagnetic.  
**Length 79 x Width 45 cm**

---

### MAIN FEATURES

Trauma	
Pediatric	
Length (cm)	79
Width (cm)	45
Weight (kg)	0,75
Usage	Transport

---

### IMAGES



---

## SEE ALSO

---



**9130**

**ME.BER. BABY FIX**

Fixing system with belts for  
transport of children inside the  
ambulance

---

---

## ACCESSORIES

---



**12034**

Red empty carrying bag for "Baby  
Fix" and "Vacuum Baby Fix"

---