



12054

Oximeb Maxi - Oxygen cylinder holder



FEATURES

Oxygen cylinders holder made of strong, antitear nylon and inside padding to prevent damages from shocks. Supplied with a carrying handle and a side pocket to contain tube and mask.

Suitable for 3 liters oxygen cylinders and, thanks to a handy fastening, it can be mounted on the side protection of the self loading stretchers and on the back rest of the sedan chairs.

MAIN FEATURES

Length (cm)	43
Width (cm)	16
Depth (cm.)	12
Weight (kg)	0,3
Material	Nylon

SEE ALSO



938
Oximeb oxygen cylinder holder

ACCESSORY OF



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF Certified self-loading stretcher



7265 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg