



7112 PROOF Evox Proof 7112

“EVOX” certified chair cot with mattress



FEATURES

The EVOX chair cot allows to place the patient in the most comfortable position during rescue and transport. Its conversion into chair allows access in all those confined spaces, such as elevators and hallways with tight corners, where a common stretcher cannot easily be manoeuvred. A gas pistons system allows to adjust EVOX in infinite positions until reaching the sitting position, relieving the operator from any sort of effort. Frame in stainless steel electro-polished, polyurethane wheels Ø 200mm, safety command to prevent accidental closure, mattress divided into four separate sections and an attractive design are the ultimate features of this device, always more required by the markets. A innovative system, consisting of four “balancing” loading wheels, allows to load Evox on different stretcher supports. Thanks to the appropriate lock in conjunction with the fixing system art. 962, the stretcher chair **EVOX art. 7112 PROOF is able to withstand shocks up to 10g in all directions, as regulated by the UNI EN 1789 norm.** The ability to accommodate patient in the most congenial position, makes Evox the ideal solution also for all the emergency waiting rooms.

MAIN FEATURES

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Height (cm)	80
Width (cm)	60
Weight (kg)	44,5
Loading capacity (kg)	220
Folded height (cm)	27
Max loading height (cm)	67
Length (cm)	193
Chair position min. length (cm)	123
Material	Stainless steel

IMAGES



ACCESSORIES



962
EN 1789 certified fixing system for
stretcher/chair "EVOX"



12100/N
Black electrowelded shaped
cushion with MeBer logo for
electrowelded mattresses



12016-E
Foldable telescopic IV pole for
stretcher EVOX Proof



608/A-MEB
2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs

STANDARD ACCESSORIES



7009/N
Black electrowelded mattress for
self loading stretcher EVOX



608/N-MEB
2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)