



## 1375

Professional bag with compartments for drugs transport "Pharma-Med"



### FEATURES

*Pharma-Med* is a new medical professional device designed by MeBer for a reliable drugs transport. This bag is made of tear-resistant, water-proof and fire-proof nylon; it is supplied with sturdy zips and it can be easily carried by a rubber coated handle.

The bottom has 4 plastic support feet to ensure maximum isolation from every kind of surface.

This device is divided into three compartments with sturdy double zips opening, allowing to separate drugs into three different categories: "*Liquid and Powder drugs*", "*Cold drugs*" and "*Blisters*".

The "*Liquid and Powder drugs*" compartment is the central one and it has three internal spaces, equipped by elastics, for holding vials with different sizes, and by two membranes made of expanded polyurethane at low memory which protect vials from breakages in case of crushing or fall.

The "*Cold drugs*", that is marked by red Cordura material, is also supplied with internal elastics for vials and reusable hot/cold storage packs (supplied with **6 Mini Dispogel**) for keeping the temperature low enough to properly store drugs.

The "*Blisters*" compartment is made of a transparent pocket with zip on one side and a transparent compartment with Velcro® opening on the other side, to allow the storage of boxes or blisters in an orderly and visible way, avoiding contamination of drugs.

*Pharma-Med* is, also, externally equipped by a useful transparent pocket for a check list sheet or a drug deadlines calendar.

Dimensions: 20 x 16 x 28 (L x H x D)

Number of vials contained: MIN 108 – MAX 120

Volume: 8960 cm<sup>3</sup>

### MAIN FEATURES

Width (cm)	20
Height (cm)	16
Depth (cm.)	28
Weight (kg)	1,1
Material	Nylon

---

IMAGES

---



---

SEE ALSO

---



**1350**

Bag with partitions for vials  
transport

---

ACCESSORIES

---



**1375-001**

Reusable heat/cold gel for drugs  
transport bags

---