



1478

Bottle for oxygen – 10 L



FEATURES

Bottles for oxygen therapy in conformity with F.U. norms, produced with stainless steel, refillable and complete of test certificate.
10 L Bottle = 200 bar - 13,25 kg. – H 955 mm – Ø 137 mm - 2000 L gaseous

MAIN FEATURES

| | |
|-------------|------|
| Height (cm) | 95,5 |
| Diameter | 13,7 |

IMAGES



SEE ALSO



1468
Bottle for oxygen – 2 L



1465
Bottle for oxygen – 0,5 L



1475
Bottle for oxygen – 7 L
