



## 7015

Folding sky mattress black for bariatric MBP platform




### FEATURES

Mattress for bariatric MBP platform made of expanded polyurethane and fireproof, waterproof and anti-tear sky cover. It is foldable into 2 parts for an easy storage and it is supplied with 6 Velcro® strips for anchoring it to the platform frame.

**Dimensions: Length x Width x Thickness: 173 x 96 x 5,5 (max thickness 10) cm**

\*All Pictures shown are for illustration purpose only.

### MAIN FEATURES

Bariatric	
Length (cm)	173
Width (cm)	96
Thickness (cm)	5,5 10

### IMAGES



---

## SEE ALSO

---



**12080-001**  
Monting kit for bariatric platform  
on self loading stretcher

---

## ACCESSORY OF

---



**12080**  
"MPB" wide bariatric platform for  
self-loading stretchers Frog  
Mercury and Gecko complete of  
belts with reel and fixing for  
stretcher



**MBS-01 MEBER BARIATRIC  
STRETCHER 7270**  
Self loading stretcher for bariatric  
transport