



## 12092

---

Battery MILWAUKEE M28 RED 5,0 Ampere for Extra Ergolift chair



---

### FEATURES

28V battery

Autonomy: 18 ramps by 20 steps with a load of 130 kg

Charging time: max 1 hour

---

### ACCESSORY OF



**Extra Ergolift 671**  
Powered stair chair



**Extra Ergolift FX 675**  
Fixed frame powered Stair Chair  
suitable for transportation



**Extra Ergolift 671**  
Powered stair chair



**Extra Ergolift FX 675**  
Fixed frame powered Stair Chair  
suitable for transportation

---