



## 7034


Fixed trolley for incubators



### FEATURES

Fixed height incubator trolley made of stainless steel structure with No. 4 wheels of  $\varnothing$  200 mm (2 with brake) for transport inside hospitals. Compact dimensions for easy manoeuvrability. This trolley is equipped with a universal interface for anchoring the incubator.

### MAIN FEATURES

Pediatric	
Material	Stainless steel
Height (cm)	74
Width (cm)	57
Length (cm)	90
Weight (kg)	30,5


---

SEE ALSO

---



**12050**  
Extended trolley incubator  
interface



**12048**  
Certified interface for incubator  
fixation for self loading trolleys

