



7300 MeBer-X

"MeBer-X" Multi-level stretcher complete with mattress and belts



FEATURES

MeBer-X, our "X" frame multilevel stretcher, has been designed in order to make the loading, unloading, lifting and lowering of the patient operations simpler and more ergonomic.

The main feature is the 9-height adjustment that facilitates transshipment and manoeuvrability in any condition. Made of aluminium alloy, the stretcher has a load capacity of up to 320 kg, which guarantees maximum safety even for the transportation of bariatric patients.

Unique features:

- 4 swivel castors with a diameter of 200 mm;
- washable and sanitisable aluminium bed;
- adjustable backrest from 0 ° to 80 °, lockable in any position;
- Trendelenburg and Fowler configuration;
- retractable loading trolley.

To be combined with the article code: art. [19162](#) MTX-2-X.

MAIN FEATURES

Weight (kg)	45,9
Length (cm)	166
Max length (cm)	197
Min height (cm)	40,5
Max height (cm)	123
Min loading height (cm)	0
Max loading height (cm)	85,5
Width (cm)	60
Loading capacity (kg)	320
Multilevel	

IMAGES



SEE ALSO



19162
MTX-2-X Double space stainless
steel stretcher support for X-frame
stretchers


ACCESSORIES



980
"MEBER-X" Multi level cot fastening system



12018N
Telescopic foldable phleboclysis stand for self loading stretchers



12120
Objects holder black for lower cart of Meber-X stretcher



12122
Objects holder black for loading trolley of Meber-X stretcher



12128
Adjustable system for headrest of Meber-X and Mercury stretchers



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



612/AN-2/MEB
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12100/N
Black electrowelded shaped cushion with MeBer logo for electrowelded mattresses



12019N
Aluminium telescopic foldable phleboclysis stand for X-frame stretcher



935
Instrument support for X-frame and Meber Max stretchers

STANDARD ACCESSORIES



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



7002

Electrowelded 4 pcs black mattress for stretchers



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)
