



675 Extra Ergolift FX

Fixed frame powered Stair Chair suitable for transportation



FEATURES

The range of Me.Ber. powered stair-chairs expands with the new Extra Ergolift FX. Thanks to the great success of the folding model Ergolift, FX borrows to grant the same high performances, but with a fixed frame, especially designed to fulfill the needs of daily transportations services or as a dedicated device in both medical and non medical facilities.

The powered chair is designed for covering the needs of modern all-round service, and it's suitable for transport on vans and ambulances. Ideal for recovering and relocating patients without efforts.

Unique features:

- wireless control panel
- telescopic handle bar
- supplied with thoracic abdominal belt
- long and high gripping tracks
- 4 dedicated points for restraint systems
- 200 mm rear wheels
- high capacity battery
- high comfort

Extra Ergolift loads up to 230 kg and can be used under rain or adverse weather conditions.

Weight 44,8 kg.

28V battery


Autonomy: 18 ramps by 20 steps with a load of 130 kg

Charging time: max 1 hour

MDD 93/42/EEC - Class I Medical Device EN 1789 Certified



MAIN FEATURES

Length (cm)	84
Max length (cm)	150
Width (cm)	55
Height (cm)	140
Weight (kg)	44,8
Loading capacity (kg)	230
Cardio	

IMAGES



SEE ALSO



Extra Ergolift 671
Powered stair chair



Extra 670/BR
Track assisted stair chair with arms



679
Narrow foldable sedan chair

ACCESSORIES



12090
Battery charger MILWAUKEE M28C
for rechargeable battery M28 RED
Extra Ergolift chair



12092
Battery MILWAUKEE M28 RED 5,0
Ampere for Extra Ergolift chair



12118
Mattress for Extra and Ergolift
chairs' transformation kit into fixed
chairs



988
Integrated control fastening
system for Extra Ergolift FX and
Proxima Pro chairs (optional)



566
Black bag with 3 belts in 2 pieces
made of black ribbon with orange
MeBer logo for Extra chairs



610
Thoracic-abdominal belt with four
fastening for padded chair



171J
Removable left lever with
ergonomic handle for Extra chairs.



172J
Removable right lever with
ergonomic handle for Extra chairs



670-004
Castor with brake diam. 200 mm
for Extra chair

STANDARD ACCESSORIES



12090
Battery charger MILWAUKEE M28C
for rechargeable battery M28 RED
Extra Ergolift chair



12092
Battery MILWAUKEE M28 RED 5,0
Ampere for Extra Ergolift chair



566
Black bag with 3 belts in 2 pieces
made of black ribbon with orange
MeBer logo for Extra chairs



12118
Mattress for Extra and Ergolift
chairs' transformation kit into fixed
chairs



610
Thoracic-abdominal belt with four
fastening for padded chair
