



12025

Black multi usage object holder for automatic loading stretcher



FEATURES

Object holder made of anti-tear nylon to be fixed under back-rest of all Me.Ber. stretchers. It is composed by a closed pocket with a zip for little objects and by a bigger open pocket for documents. It is fixed to the trolley tube of the stretcher with practical Velcro® belts. **Dimensions: cm 50x37**

MAIN FEATURES

Width (cm)	50
Height (cm)	37

ACCESSORY OF



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with variable height and four swivel castors



Mercury 7070/4RG PROOF
Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury Lite 7082/4RG Proof
Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury 7072/4RG PROOF
Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



Gecko 7122/AZ
"GECKO" self loading stretcher entry-level light blue



Gecko 7122/V
"GECKO" self loading stretcher entry-level green



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium stretcher with trendelenburg certified



Frog 7210
"FROG" yellow self loading stretcher with trendelenburg



Frog 7210 PROOF
"FROG" yellow self loading stretcher with trendelenburg certified



Frog 7240 PROOF
"FROG PLUS" self loading stretcher with trendelenburg certified



MBS-01 MEBER BARIATRIC STRETCHER 7270
Self loading stretcher for bariatric transport



Winner 910
Self loading stretcher



Winner 910 PROOF
Self loading stretcher - certified



Winner 914
Self loading stretcher with trendelenburg and fowler



Winner 914 PROOF
Certified self loading stretcher with Trendelenburg and Fowler



Mercury Lite Cinque 7095/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version
