



11220 Oxytest

OXYTEST portable pulse oximeter complete with adult sensor, USB cable and AAA battery



FEATURES

Oxytest is a portable pulse oximeter, easy to use and light for transport. Complete with:

- Oled display, four directions to display arterial oxygen saturation, pulse rate and pulse strength
- battery voltage indicator
- standby or automatic suspension

Thanks to the sensors art. 11224 and art. 11222 (not included in the package), it is possible to use Oxytest also for children and new born.

MAIN FEATURES

Height (cm)	14,5
Width (cm)	8,3
Depth (cm.)	3,3
Weight (kg)	0,24
Cardio	

IMAGES



SEE ALSO



16500 Sterilair Pro Biological air treatment device -Wall version



16502 Sterilair Pro Biological air treatment device -Stand version



Prevention and individual protection



8377 Infrared thermometer



8506Pulse oximeter Oxy-110



11200 Vital Test Handheld pulse oxymeter complete with adult sensor

ACCESSORIES



11224 Pediatric sensor for Oxytest



11222 New born sensor for Oxytest

STANDARD ACCESSORIES



11226Adult sensor for pulse oximeter
Oxytest