



8432

Teca di colore bianco mod. standard per defibrillatore



FEATURES

Metal cabinet to contain a Mobile Emergency Defibrillator. Enamelling with polyester powder for external RAL9010 cooked at 180° One door with transparent inspection window, with an easy zinc-coated S/S sealed lock, anti dust roof, back hanging wall holes.

SEE ALSO



8396

Defibrillator metal cabinet with thermoregulation and alarm for outer use

ACCESSORIES



8432-001
Pedestal for art. 8432

ACCESSORY OF



8440
Semiautomatic defibrillator mod.
Technoheart Plus AED



8442
Semiautomatic defibrillator mod.
Technoheart Plus AED Trainer
