



707

Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



FEATURES

Slash is MeBer's innovative mattress protector transfer sheet. Excellent ally for moving the patient from bed to stretcher and vice versa. Equipped with 6 handles and made of fireproof, tearproof and waterproof nylon.

Placed on self-loading stretchers, it easily lends itself to immediate sanitation between one patient and another. The reduced sizes are specially designed to facilitate operations during the handling phases.

Compatible with all stretcher models.

SLASH IS NOT A CARRYING SHEET.

MAIN FEATURES

| Length (cm) | 188 |
|-----------------------|-----|
| Width (cm) | 50 |
| Weight (kg) | 1 |
| Loading capacity (kg) | 160 |

IMAGES



ACCESSORY OF



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels self-

Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury Lite 7080/4RG Proof Ambulance cot in aluminium with variable height and four swivel



Mercury 7070/4RG PROOF Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version



7265 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg