



## 12084

---

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



---

### FEATURES

The bariatric wings are a practical accessory that allows you to convert the self-loading stretcher into a practical bariatric recovery tool. Thanks to the small footprint (and a weight of only 14 kg) they can be easily transported by ambulance to be used when needed. Their use allows to face in safety, and in the most adequate way possible, every recovery scenario in the presence of bariatric patients. The kit includes two wings that extend the support surface of the stretcher by over 40 cm; assembly and disassembly on the self-loading stretcher are extremely simple and intuitive and do not require the use of tools. If mounted on our MBS01 they allow an effective load capacity of 335 kg. The bariatric wings are equipped with mattresses (art. [12170](#) and art. [12172](#)) and belts with reel. They can be purchased in kits including both the right and the left (art. 12084) or separately (art. 12086 Right bariatric wing and 12088 left bariatric wing).

**IMPORTANT:** For Mercury self loading stretcher use art. [12200](#).

Dimensions:  
1365x205x85 mm (a wing)

---

### MAIN FEATURES

Length (cm)	136,5
Width (cm)	20,5
Depth (cm.)	8,5

---

## IMAGES

---



---

## SEE ALSO

---



**12200**

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

---

## ACCESSORIES

---



**12084-001**

Cover for bariatric wings set



**12170**

Black mattress for right bariatric wing for self-loading stretchers



**12172**

Black mattress for left bariatric wing for self-loading stretchers



**12080-001**

Monting kit for bariatric platform on self loading stretcher

---

---

## ACCESSORY OF

---



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg



**Frog 7210**  
"FROG" yellow self loading stretcher  
with trendelenburg



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified



**MBS-01 MEBER BARIATRIC  
STRETCHER 7270**  
Self loading stretcher for bariatric  
transport



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium  
stretcher with trendelenburg  
certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher  
with trendelenburg certified

---