



19136-956

"MX-956" fixed single space stretcher support with ramp of mm. 180, equipped with en1789 certified fixing system art. 956 - EN1789 certified



FEATURES

Fixed single space stretcher support with ramp mm 180. On request supplied with UNI EN 1789 certified blocking devices Art. 956.



MAIN FEATURES

| | |
|-----------------------|-------|
| Width (cm) | 60,5 |
| Length (cm) | 208,5 |
| Weight (kg) | 50 |
| Loading capacity (kg) | 400 |
| Height (cm) | 13,5 |

ACCESSORIES



19114
130 mm ramp for stretcher supports



19112
100 mm ramp for stretcher supports