




# 19146-955

"MX-2-955" double space fixed stretcher support with ramp of mm. 350, equipped with certified fixing system art. 955 - EN1789 certified



## FEATURES

Double space fixed stretcher support with ramp mm 350. On request supplied with UNI EN 1789 certified blocking devices Art. 955.  
Gas spring included.



## MAIN FEATURES

Width (cm)	63
Length (cm)	223
Weight (kg)	68,2
Loading capacity (kg)	400
Height (cm)	21,5

---

ACCESSORIES

---



**19110**

Stainless steel plate mm 1940 x 572  
for stretcher supports



**787L**

Gas spring for ramp

---

STANDARD ACCESSORIES

---



**821G**

Gas spring for ramp