




## 19146-955

---

"MX-2-955" double space fixed stretcher support with ramp of mm. 350, equipped with certified fixing system art. 955 - EN1789 certified



---

### FEATURES

Double space fixed stretcher support with ramp mm 350. On request supplied with UNI EN 1789 certified blocking devices Art. 955.  
Gas spring included.



---

### MAIN FEATURES

---

Width (cm)	63
Length (cm)	223
Weight (kg)	67,4
Loading capacity (kg)	400
Height (cm)	21,5

---

## ACCESSORIES

---



**19110**

Stainless steel plate mm 1940 x 572  
for stretcher supports



**787L**

Gas spring for ramp

---

---

## STANDARD ACCESSORIES

---



**821G**

Gas spring for ramp

---