




6100


Long folding stretcher



FEATURES

Stretcher that folds in two, equipped with feet and handles for transport. Aluminium frame and PVC impermeable and tear-resistant bed. Ideal for in-surgery transport, easily positioned. Available in two different dimensions.

MAIN FEATURES

Foldable	
Lying surface	Without Trendelenburg
Length (cm)	167
Width (cm)	50
Height (cm)	13
Folded height (cm)	83
Thickness (cm)	13
Weight (kg)	5,5

SEE ALSO



6300/B
Vacuum mattress for veterinary use



Maxima 6200
Pediatric anti trauma scoop
stretcher
