



614

Certified horizontal fastening system for scoop stretchers



FEATURES

Aluminium support with safe quick-locking device for wall mount or horizontal storage of a scoop stretcher.



Device certified according to the European harmonized standards UNI EN 1789

MAIN FEATURES

Length (cm)	50,5
Width (cm)	2,5
Height (cm)	6,5
Weight (kg)	0,7


SEE ALSO



Maxima 630/A
Orange painted anti trauma scoop
stretcher



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 630
Grey anodized anti trauma scoop
stretcher