



630/A Maxima

Orange painted anti trauma scoop stretcher



FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

MAIN FEATURES

| | |
|-----------------------|---|
| Size | Adults |
| Usage | Rescue |
| Height (cm) | 7 |
| Width (cm) | 42 |
| Length (cm) | 167 |
| Max length (cm) | 202 |
| Folded length (cm) | 120 |
| Epaisseur pliée (cm) | 9 |
| Weight (kg) | 10 |
| Loading capacity (kg) | 170 |
| 10G Certification |  |
| Material | Painted aluminium |

SEE ALSO



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630/N
Black painted anti trauma scoop
stretcher



Maxima 630
Grey anodized anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher

ACCESSORIES



616
Certified vertical fastening system
for scoop stretchers



616/A
Certified vertical fixing system for
scoop stretchers



614
Certified horizontal fastening
system for scoop stretchers



630/CS
Carrying case for scoop stretchers



Double 624
Head immobilizer for scoop
stretcher



603/MEB
Nylon bag containing 3 belts Art.
607/Ameb