



## 660 Folding stretcher

Collapsible stretcher foldable in 4 pieces, telescopic handles



### FEATURES

This 4 pieces foldable stretcher is the most complete and versatile in the range. The sturdy tubular aluminium frame makes it light and robust. The stretcher features 4 telescopic handles that allow for a larger support surface. The canvas is in fire-resistant impermeable, tear-resistant, easily cleanable grey nylon. Supplied with 2 quick-release straps for patient immobilization (Art. 605/Nmeb).

### MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Lying surface	Without Trendelenburg
Min height (cm)	18
Epaisseur pliée (cm)	21
Width (cm)	57
Length (cm)	192
Max length (cm)	216
Folded length (cm)	94
Weight (kg)	9,3
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

### IMAGES



---

## SEE ALSO

---



**Folding stretcher 1550**  
Military folding stretcher with  
green canvas

---

---

## ACCESSORIES

---



**660/CS**  
Carrying case for collapsible  
stretcher Art. 660

---

---

## STANDARD ACCESSORIES

---



**605/Nmeb**  
50 mm black belt in one piece with  
plastic buckle and ME.BER. ribbon

---