



## 7014/A

---

4 pcs. Skay fire-proof electrowelded mattress



---

### FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, **can be used on almost all MeBer self loading stretchers (not suitable for stretcher Evox).**

Easy to clean and sanitize.

**Dimensions: cm 186 x 47 x 9.**

\*All Pictures shown are for illustration purpose only

---

### MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

---

## ACCESSORY OF

---



**Fixed stretcher 1515 PROOF**  
Fixed stretcher with rigid bed



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels -  
high version



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels