



## 8004

---

Colored aneroid sphygmomanometer with carrying case



---

### FEATURES

Range 0-300 mm Hg. Accuracy +/-3 mm Hg. Manometer NON STOP PIN. Adult Velcro® cuff and metal ring. Standard model in vinyl zipper case.

---

### MAIN FEATURES

Size

Adults

---

### SEE ALSO



**1740**  
Rotating wall aneroid  
sphygmomanometer



**8012**  
Red palm-style aneroid  
sphygmomanometer plus

---

## ACCESSORIES

---



**8143**  
Double tube pediatric cuff



**8142**  
Single tube pediatric cuff



**8140**  
Adults cuff with single tube



**8141**  
Adults cuff with double tube

---