



8004

Colored aneroid sphygmomanometer with carrying case



FEATURES

Sphygmomanometers with nylon cuff with velcro closure, reliable manometer with clip to be fixed on the cuff. Scale 0-300 mmHg. From 260 mmHg to 15 mmHg max deflation time 10 sec. Chromium-plated air tap with air release regulation by screw valve. Supplied with a simulated-leather soft pouch with zip. Latex free.

Range 0-300 mm Hg. Accuracy +/-3 mm Hg. Manometer NON STOP PIN. Adult Velcro® cuff and metal ring. Standard model in vinyl zipper case.

MAIN FEATURES

Size

Adults

SEE ALSO



1740
Rotating wall aneroid
sphygmomanometer



8012
Red palm-style aneroid
sphygmomanometer plus

ACCESSORIES



8143
Double tube pediatric cuff



8142
Single tube pediatric cuff



8140
Adults cuff with single tube



8141
Adults cuff with double tube
