



8006

Aneroid sphygmomanometer



FEATURES

Sphygmomanometers with nylon cuff with velcro closure, reliable manometer with clip to be fixed on the cuff. Scale 0-300 mmHg. Floating zero, nylon velcro cuff, latex bulb with chromium plated valve. Chromium-plated air tap with air release regulation by screw valve. Supplied with a simulated-leather soft pouch with zip.

MAIN FEATURES

Size

SEE ALSO





8012Red palm-style aneroid sphygmomanometer plus



1740 Rotating wall aneroid sphygmomanometer

ACCESSORIES









ACCESSORY OF



1366Red knapsack "115" with side pockets and compartments COMPLETE