



## 8009

---

Shock proof aneroid sphygmomanometer for children



---

### FEATURES

Manometer range 0-300 mm/Hg. Accuracy +/- 3 mm/Hg. Paediatric Velcro® cuff and metal ring. Standard model in vinyl with closure zipper case.

---

### MAIN FEATURES

---

Size

Children

Pediatric



---

SEE ALSO

---



**8012**  
Red palm-style aneroid  
sphygmomanometer plus



**8008**  
One hand aneroid  
sphygmomanometer



**8004**  
Colored aneroid  
sphygmomanometer with carrying  
case



**8006**  
Aneroid sphygmomanometer



**8011**  
Set newborn/pediatric aneroid  
palmar sphygmomanometer

---