



## 8020/N

---

Black flat phonendoscope



---

### FEATURES

Flat phonendoscope with anodized stamped and lathed aluminium head diameter 46 mm. Threaded metal ring for membrane blocking. Metal chromium plated binaural, length 150 mm. Y rubber tube in one piece, diameter 5 mm. length 60 cm. Single printed carton box. Available in 3 different colors.

---

### MAIN FEATURES

Size	Adults
Cardio	

---

## SEE ALSO

---



**8025/N**  
Black stethophonendoscope



**8025/G**  
Yellow stethophonendoscope



**8025/R**  
Red stethophonendoscope



**8020/R**  
Red flat phonendoscope



**8020/G**  
Yellow flat phonendoscope



**8030/B**  
Blue stethophonendoscope  
"Rappaport"

---

## ACCESSORY OF

---



**1366**  
Red knapsack "115" with side  
pockets and compartments  
COMPLETE

---