



862 Sirius

Rigid cervical collar in 1 piece size PEDIATRIC



FEATURES

Rigid cervical collar in 1 piece made with transparent polyethylene with E.V.A. washable inner lining fixed by plastic rivets. Adjustable rip-off fastening which different colours indicate the various sizes. Central hole for tracheotomy. Available in 6 different sizes.

MAIN FEATURES

| | |
|-----------|----------------|
| Size | Children |
| Usage | Immobilization |
| Trauma | |
| Pediatric | |

SEE ALSO



Antares 854

Cervical collar set 4 pieces with carrying case Art. 842 - 844 - 846 - 848



Sirius 876

Cervical collar set 6 pieces with carrying case



Antares 856

Cervical collar set 6 pieces with carrying case



Sirius 874

Cervical collar set 4 pieces with carrying case Art. 862 - 864 - 866 - 868

ACCESSORY OF



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE
