



## 880 Halley

Vacuum splint for leg



### FEATURES

Produced with absolute reliability materials this splints allow the perfect immobilization of the lower and upper limbs by taking away the inside air. They are supplied with suction valve positioned on the out-side to avoid problems to the patient. It is possible to order the single vacuum splint or the complete set including splints for leg, arm, forearm, little pump for the vacuum and nylon bag. All splints are X-ray transparent.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	109
Width (cm)	67
Material	PVC

### IMAGES



---

## SEE ALSO

---



**Halley 886**  
Vacuum splint set with pump and carrying case



**Halley 884**  
Vacuum splint for forearm



**Halley 886**  
Vacuum splint set with pump and carrying case

---

## ACCESSORIES

---



**889**  
Little aluminium pump