



## 9150 X-Lock

Ergonomic head immobilizer for Ergon stretcher



### FEATURES

One piece head immobilizer with integrated head strap, for a fast head immobilization directly on the stretcher. Made with laminated polyesters and bi-density polyurethane foams.

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	43
Width (cm)	27
Height (cm)	140

### SEE ALSO



**Ergon 16150/V**  
Variable geometry comfort  
stretcher - green/blue



**Ergon 16150/G**  
Variable geometry comfort  
stretcher - yellow/blue

---

## ACCESSORIES

---



**Art. 9150-001**  
Blue carrying case for X-Lock  
ergonomic head immobilizer

---