



938

Oximeb oxygen cylinder holder



FEATURES

Cylinder holder produced with strong, antitear nylon and inside padding to prevent damages from knocks. It is supplied with a handle for the transport and a side pocket for tube and mask.

Can be used with cylinders of 2 liters and, thanks to a handy fastening, it can be mounted on the side protection of the automatic loading stretchers and on the back rest of the sedan chairs.

MAIN FEATURES

| Weight (kg) | 0,3 |
|-------------|-------|
| Length (cm) | 31 |
| Width (cm) | 19 |
| Depth (cm.) | 13 |
| Material | Nylon |

SEE ALSO



12054Oximeb Maxi - Oxygen cylinder
holder

ACCESSORY OF



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF Certified self-loading stretcher



7265 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg