



939

Instrument support for stretchers



FEATURES

Instrument support to be fixed on Me.Ber. stretcher frames' to transport medical equipments, defribillators and monitors. It is supplied with belts to fix devices and is foldable to be stored inside ambulances. Dimensions: open cm 34x56x37h, closed hight cm 8

MAIN FEATURES

Width (cm)	56
Depth (cm.)	34
Height (cm)	37
Epaisseur pliée (cm)	8
Weight (kg)	5,6
Loading capacity (kg)	30

IMAGES





SEE ALSO



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels

ACCESSORY OF



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury EL 7094 PROOF Certified self-loading stretcher high version



7265 Proof Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg