



## 950/2

---

Rear fixing system for roll-in H=160 mm



---

### FEATURES

Built in stainless steel fusion. It is mounted on the ambulance platform to safely fix the automatic loading stretcher. It is available in 4 different heights.

---

### MAIN FEATURES

Height (cm)

16

---

### SEE ALSO



#### 950/1

Rear fixing system for roll-in H=188 mm



#### Rear fixing system

Rear fixing system for roll-in H=148 mm



#### Rear fixing system

Rear fastening system for ambulance cots H=129 mm

---