



955 PROOF

EN 1789 fixing system for Proof roll-in stretchers



FEATURES

Designed to guarantee maximum stability of the stretcher on board the rescue vehicle, they are made entirely of stainless steel for longer life and corrosion resistance. Their easy-to-install structure ensures that the stretcher or trolley, and patient are securely restrained in all directions. Art. 955 Proof is suitable for fixing all our 10G. roll-in stretchers whilst art. 956 Proof is suitable for fixing "Titanus" trolleys.



Device certified according to the European harmonized standards UNI EN 1789.

| Height | Width | Lenght | Weight |
|--------|-------|--------|--------|
| 13 cm | 48 cm | 177 cm | 7,7 kg |

MAIN FEATURES

| Height (cm) | 13 | |
|-------------|-----|--|
| Length (cm) | 177 | |
| Weight (kg) | 7.7 | |

SEE ALSO



956 PROOF
Fixing system for Titanus automatic loading trolley

ACCESSORIES



955-001 Front reverse-plate for 955-956 Proof certified fixing system



955-002 Rear reverse-plate for 955-956 Proof certified fixing system

ACCESSORY OF



MBS-01 MEBER BARIATRIC STRETCHER 7270 Self loading stretcher for bariatric transport



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading

stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF

Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite 7080/4RG Proof

Ambulance cot in aluminium with variable height and four swivel castors



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury 7070/4RG PROOF

Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF

Certified stainless steel variable height stretcher high version with 4 rotating wheels



Frog 7240 PROOF "FROG PLUS" self loading stretcher with trendelenburg certified



Frog Lite 7260 PROOF "FROG LITE" self loading aluminium stretcher with trendelenburg certified



Frog 7210 PROOF "FROG" yellow self loading stretcher with trendelenburg certified



Mercury Lite Cinque 7095/4RG

ProofCertified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG

ProofCertified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version



Mercury Lite 7082/4RG Proof

Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury EL 7094 PROOF Certified self-loading stretcher

high version



7265 Proof

Aluminium self loading stretcher "Frog lite Advance" with exctractable handles and Trendelenburg