



## 16160/G Ergon Junior

Pediatric variable geometry comfort stretcher - yellow/blue



### FEATURES

[Visit Ergon website](#)

The greatest product for the smallest patients

ERGON shares its advanced technology with the junior model. All benefits and main features of the adult stretcher are preserved in its Junior version.

The **Flexilock locking system** was moved from the centre to the side, in order to allow operations without invading the head space. The large blade allows for maximum comfort even to pediatric patients.



### MAIN FEATURES

Size	Children
Usage	Rescue Immobilization Transport
Height (cm)	7,5
Width (cm)	43
Length (cm)	129
Max length (cm)	153
Folded length (cm)	85
Epaisseur pliée (cm)	9
Weight (kg)	8,4
Loading capacity (kg)	170
Type	
Pediatric	<input checked="" type="checkbox"/>
Trauma	<input checked="" type="checkbox"/>
Alpine	<input checked="" type="checkbox"/>
Material	HC3 - High Comfort Composite Cell Aluminium alloy

---

## SEE ALSO

---




**Ergon Junior 16160/V**  
Pediatric variable geometry  
comfort stretcher - green/blue

---

---

## ACCESSORIES

---



**X-FIX 9156**  
Restraint system with straps and  
carabiners



**959**  
Wall mount system for variable  
geometry stretcher Ergon



**16150-002**  
Storage bag for Ergon stretcher



**16150-001**  
Winch bag for Ergon

---

---

## STANDARD ACCESSORIES

---



**584**  
Black bag with 3 orange straps in 2  
pcs and quick-release plastic  
buckle for Ergon stretcher

---