



7013/A

3 pcs. Skay fire-proof electrowelded mattress



FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, **can be used on almost all MeBer self loading stretchers with Trendelenburg including Cronos (not suitable for stretcher Evox and stretcher 1515).**

Dimensions: cm 186 x 47 x 9.

*All Pictures shown are for illustration purpose only

MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

SEE ALSO



7013/N
3 pcs. Skay fire-proof electrowelded
mattress

ACCESSORY OF



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury 7072/4RG PROOF
Certified stainless steel variable
height stretcher high version with 4
rotating wheels



Winner 910 PROOF
Self loading stretcher - certified



Winner 914 PROOF
Certified self loading stretcher with
Trendelenburg and Fowler



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels



Mercury EL 7092 PROOF
Certified self-loading stretcher
