



## 430

Arcos rigid splint set with carrying bag



### FEATURES

New rigid splints, made with fire-proof and water proof PVC, allow a correct immobilization of the limbs. They are X-ray translucent. Equipped with 4 Velcro® closures with adjustable length allowing operator to apply them easily and safely. Shaped to strongly and easily fix to ankle and upper part of the leg, they are supplied with foot support integrated.

Available in two sizes:

Adult arm/ pediatric leg: length cm. 52 - base cm 9

Adult leg: length cm 74 - base cm 13

### MAIN FEATURES

Trauma



Usage

Immobilization

### IMAGES



---

## SEE ALSO

---



**427**  
Arcos short rigid splint for upper  
limbs



**425**  
Arcos long rigid splint for lower  
limbs

---