




629 Jelly


Orange standard head immobilizer for spinal board



FEATURES

Adjustable to every models of spinal boards. It is made of expanded materials at closed cells covered by a vinyl material.

MAIN FEATURES

| | |
|-------------|---|
| Trauma |  |
| Width (cm) | 40,5 |
| Depth (cm.) | 25,5 |
| Height (cm) | 16,5 |

ACCESSORY OF




Galaxy 9010
"GALAXY" Yellow spinal board



Galaxy 9012
"GALAXY" Orange spinal board



Galaxy 9016
"GALAXY" Red spinal board



Galaxy-P 9040
Yellow spinal board with pins



Galaxy-P 9042
Orange spinal board with pins



Niliix 9050
Yellow basic spinal board with pins