




629 Jelly


Orange standard head immobilizer for spinal board



FEATURES

Adjustable to every models of spinal boards. It is made of expanded materials at closed cells covered by a vinyl material.

MAIN FEATURES

Trauma	
Width (cm)	40,5
Depth (cm)	25,5
Height (cm)	16,5

ACCESSORY OF




Galaxy 9010
"GALAXY" Yellow spinal board



Galaxy 9012
"GALAXY" Orange spinal board



Galaxy 9016
"GALAXY" Red spinal board



Galaxy-P 9040
Yellow spinal board with pins



Galaxy-P 9042
Orange spinal board with pins



Nilix 9050
Yellow basic spinal board with pins