



## 15244

Folding Scangliss



### FEATURES

Foldable transfer board to use under the patient allowing a safe transfer from stretcher to bed. It is translucent, a-magnetic and it is divided in two sections. This device consists of a rigid polyethylene core, padded with an amorphous material of 1,5 cm thickness, all wrapped in a sliding nylon sheet. Scan Gliss is CE certified according to Directive 93/42 EEC Medical Devices.

**Maximum load allowed: 180 kg.**

### MAIN FEATURES

|                       |              |
|-----------------------|--------------|
| Length (cm)           | 180          |
| Width (cm)            | 50           |
| Height (cm)           | 2            |
| Weight (kg)           | 4            |
| Loading capacity (kg) | 180          |
| Material              | Polyethylene |

### SEE ALSO



**1398**

White/blue absorbing and transfer sheet cm. 220x100



**15180**

Big Gurney

---

## ACCESSORIES

---



**15244-001**  
Nylon cover for Scangliss

---