



15244

Folding Scangliss



FEATURES

Foldable transfer board to use under the patient allowing a safe transfer from stretcher to bed. It is translucent, a-magnetic and it is divided in two sections. This device consists of a rigid polyethylene core, padded with an amorphous material of 1,5 cm thickness, all wrapped in a sliding nylon sheet. Scan Gliss is CE certified according to Directive 93/42 EEC Medical Devices.
Maximum load allowed: 180 kg.

MAIN FEATURES

Length (cm)	180
Width (cm)	50
Height (cm)	2
Weight (kg)	4
Loading capacity (kg)	180
Material	Polyethylene

SEE ALSO



15180
Big Gurney

ACCESSORIES



15244-001
Nylon cover for Scangliss
